BIOLÁTEX ASO

REF K004

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for the qualitative and semi-quantitative determination of Antistreptolysin O by the agglutination of latex particles without dilution of the sample. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Latex

This method is based on a particle agglutination reaction of latex coated with antibody, especially treated to prevent non specific agglutination. Agglutination is visible in samples with concentrations equal to or greater than 200 UI/mL, according to International Standards established by the references of the WHO.

REAGENTS

1 - Number 1 - ASO Latex - Store between 2 and 8°C. **Do not freeze.** Contains: Latex Particles in suspension sensitized with Streptolysin O.

2 - Number 2 - Positive Control - Store between 2 and 8°C. Contains: Latex Agglutinating Solution.

3 - Number 3 - Negative Control - Store between 2 and 8°C. Contains: 0,85% Saline Solution.

PRESENTATION

Reagent	Presentation	
	1 (Normal Packing)	2 (Economic Packing)
Reagent N° 1	2 mL	2 mL
Reagent N° 2	1 mL	--
Reagent N° 3	1 mL	--

The normal package is accompanied by spatulas and test cards.

EQUIPMENTS AND OPERATIONAL INPUTS

Slides or dark bottom plates, spatulas, watch or stopwatch. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Do not freeze reagents.
- 5- Always use reagents from same lot.

6- Do not use lipemic serum. Do not use plasma.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum, with no previous dilution. The analyte is stable for 02 days if kept between 2 and 8°C.

PROCESS DESCRIPTION

TECHNIQUE

In each circle of the slide place the following:

	Circle N° 1	Circle N° 2	Circle N° 3
Negative Control	20 µL	--	--
Positive Control	--	20 µL	--
Sample	--	--	20 µL
Reagent N° 1	20 µL	20 µL	20 µL

Homogenize with the aid of a spatula using the entire length of each circle of the slide. Soon after, agitate the slide with circular motions for two minutes. Make the reading with an artificial light, using a dark background to facilitate test interpretation.

A clear agglutination indicates the presence of Antistreptolysin O at a concentration equal to or greater than 200 UI/mL. In this case, perform the semiquantitative test.

SEMIQUANTITATIVE TEST

1- Perform sample dilutions with saline, starting from initial sample (1:2, 1:4, 1:8, 1:16, 1:32, etc.).

2- Follow the process as described in qualitative test for each of the dilutions.

It will be considered as title the biggest serum dilution that presents agglutination.

RESULTS

Positive: Clear presence of agglutination.

Negative: Absence of agglutination (homogenized suspension).

CALCULATIONS

Sample	Concentration (UI/mL)
No dilution	200
1 / 2	400
1 / 3	600
1 / 4	800
1 / 5	1000

Results can be expressed in titles or in UI/mL.

UI/mL = 200 x title from last dilution (dilution number)

Negative Tests: Express the results with negative or minor than 200 UI/mL.

PROCEDURE LIMITATIONS

Do not use plasma, hemolysed or lipemic serum, as they may produce nonspecific agglutination. By correlating methods for determining Anti-streptolysin O, to verify the sensitivity of the reagents.

The results should only be compared when expressed in UI/mL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's sensitivity was determined through the reference material from WHO.

REFERENCE VALUES

Up to 200 UI/mL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Biolatex ASO kit was compared with another method to measure the presence of Antistreptolysin O.

According to results from 10 clinical samples, all methods showed excellent correlation. With these results we can conclude the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 5 samples with different concentrations, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Results Obtained
01	20	800	800
02	20	200	200
03	20	Negative	Negative
04	20	Negative	Negative
05	20	400	400

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 samples with different concentrations, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Results Obtained
01	20	Negative	Negative
02	20	Negative	Negative
03	20	800	800
04	20	400	400
05	20	Negative	Negative

Sensitivity

Clinical or diagnostic sensitivity and methodological sensitivity were determined through samples of 20 human serum truly positive. Tests show 100% positive results, no false negative results were found.

DIAGNOSTIC SIGNIFICANCE

The Streptolysin O is a toxic immunogenic exoenzymes produced by many β -hemolytic streptococci group A. Taking into account the large amount of exoenzymes released *in vivo*, the determination the antibody response against the presence of Streptolysin O has become a routine procedure for diagnosis and treatment of rheumatic fever, acute glomerulonephritis, scarlet fever, tonsillitis, erysipelas, puerperal sepsis and other infections streptococcus group A.

NUMBER OF TESTS

Presentation 1 (Normal Packing) 100 tests / 20 μ L of Sample / 20 μ L of Reagent

Presentation 2 (Economic Packing) 100 tests/20 μ L of Sample / 20 μ L of Reagent

BIBLIOGRAPHIC REFERENCES

- 1 - PLOTZ and SINGER, A. J. Med. 22, 1979.
- 2 - ADAMS, L. E.; HESS, E., J. Amer. Technol, 48, 1978.
- 3 - NORMAUSSELL, D., Immunochemistry, 9, 1972.
- 4 - DITO, W., Am. Soc. Clin. Pat., 69, 1976.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Biolatex ASO kit: 10269360137

Review: March/2013

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED